



Corneal OCT for Refractive & Cataract Surgeons

David Huang, MD, PhD
 Peterson Professor of Ophthalmology
 Professor of Biomedical Engineering
 Casey Eye Institute, Oregon Health & Science University
 Portland, Oregon

CASEY EYE Institute
OHSU

Financial Interests:
 OHSU and Dr. D. Huang have a significant financial interest in Optovue, a company that may have a commercial interest in the results of this research and technology. These potential conflicts of interest have been reviewed and managed by OHSU. Optovue, Inc.: patent royalty, equipment loan, stock ownership.


1

Corneal OCT Applications

CORNEAL PACHYMETRY AND EPITHELIUM MAPPING

2

Pachymetry Mapping with OCT



Li Y, et al. (Huang D) *Ophthalmology* 2006;113(5):792-9

3

OCT Pachymetry Agreement with Other Modalities

N	OCT	OCT CCT (μm)	OCT- Other modalities (μm)		
			Scheimpflug (Pentacam [®] or Galle [®])	Slit-scanning (Orbscan II)	Ultrasound (Sonogage [®] or Sonomed [®])
50	RTVue ¹	536.9	--	-0.3±12.1	-19.7±10.5 ^C
66	RTVue ²	532.8	-6.0±4.8 ^A	--	--
50	Casia ³	547.2	-11.7±6.0 ^B	-7.2	-9.2 ^D

CCT = central corneal thickness

1. Li Y, et al. (Huang D) *J Cataract Refract Surg* 2010;36(5):826-831.
 2. Huang J, et al. (Wang Q) *PLoS One* 2014;9(5):e98316.
 3. Lee YW, et al. (Choi CY) *J Cataract Refract Surg* 2015;41(5):1018-1029.

4


OCT Pachymetry Repeatability and Reproducibility

	OCT	Central D<2mm	D=2~5mm
Repeatability	RTVue ¹	1.3 μm	1.8~3.8 μm
	RTVue ²	2.1 μm	2.9~5.5 μm
	Avanti ³	1.3 μm	3.1~6.7 μm
Reproducibility	Casia ⁴	2.5 μm	3.8~6.1 μm
	RTVue ⁵	2.1 μm	3.6 μm

1. Li Y, et al. (Huang D) *J Cataract Refract Surg* 2010;36(5):826-831.
 2. Huang J, et al. (Wang Q) *Ophthalmology* 2013;120(10):1951-1958.
 3. Unpublished data.
 4. Neri A, et al. (Neri A) *Acta Ophthalmol* 2012; 90:e452-e457.
 4. Prakash G, et al. (Agarwal A). *Am J Ophthalmol* 2009;148(2):282-290 e282.

5

Pachymetry Mapping with OCT



Optovue Avanti

Avanti OCT 6-mm Pachymetry Map
 "Pachymetry+CPwr" scan pattern

Pachymetry

SN-I(2-5mm): [22] S-I(2-5mm): [23]
 Min: [522] Location Y: [-159]
 Min-Median: [-23] Min-Max: [-72]
 Min thickness at(-0.312mm, -0.159mm) indicated as *

Normal cornea is:
 - Mostly green on the pachymetry map
 - Thinnest cornea locates inside the central D=2 mm circle

David Huang, MD, PhD www.COOLab.net

6

OCT Keratoconus Pachymetric Parameters

Parameter	Explanation
IT-SN	Average thickness of the IT octant minus that of the SN octant
I-S	Average thickness of the inferior (I) octant minus that of the superior (S) octant
Min	Minimum corneal thickness
Min - Med	Minimum corneal thickness - median corneal thickness
Y Min	Y coordinate of minimum corneal thickness

1. Li Y, et al. (Huang D) *Ophthalmology* 2008;115(12):2159-2166.
 2. Qin B, et al. (Huang D) *J Cataract Refract Surg* 2013;39(12):1864-1871.

7

OCT Pachymetry Based Keratoconus Risk Score

Keratoconus Risk Score Table
 Available for downloading @ <http://www.coollab.net/resources>

Patient Name: _____

Variables (µm)	0	1	2	3	OD	OS
SN-IT	<33	33-42	43-51	>51		
Minimum	>499	499-476	475-455	<455		
Minimum-Median	>-21	-21--25	-26--29	<-29		
S-I	<30	30-40	41-49	>49		
Ymin	>734	-734--1069	-1070--1353	<-1353		
Keratoconus Risk Score						

Keratoconus risk:
 Keratoconus risk score 0-3: low risk, ≥4: high risk.

- Each variable will be assigned a score of 1, 2, 3 if it exceeds 20, 5, 1 percentile thresholds.
- The keratoconus risk score of the eye is the summation of all scores.

Qin B, et al. (Huang D). *J Cataract Refract Surg* 2013;39(12):1864-71.

8

Corneal Epithelial Imaging with FD-OCT

Li Y, et al. (Huang D) *Ophthalmology* 2012; 119(12):2425-33.

9

OCT Epithelial Thickness Map

iVue OCT 6-mm Epithelial Thickness Map
 "Cornea Pachymetry" scan pattern

Epithelium	
Epithelium statistics within central 5 mm	
S (2-5mm):	51 I (2-5mm): 53
Min:	50 Max: 55
Std Dev:	1.4 Min-Max: -5
Min/Max thickness indicated as +/-	

Normal cornea is:

- Mostly green on the epithelial map
- Superior slightly thinner than the inferior
- Epithelial Std Dev < 3.6

David Huang, MD, PhD www.COOL-Lab.net

10

OCT Epithelial Map Repeatability

OCT Device	Central D<2mm	D=2~5mm
RTVue ¹	0.7 µm	0.7~1.1 µm
RTVue ²	0.7 µm	0.6~0.9 µm
Avanti ³	1.6 µm	1.2~1.7 µm
Cirrus HD ⁴	1.5 µm	1.3~1.5 µm


1. Li Y, et al. (Huang D) *Ophthalmology* 2012;119(12):2425-2433.
 2. Ma XJ, et al. (Koch DD) *Cornea* 2013;32(12):1544-1548.
 3. Hashmani N, et al. (Hashmani S) *Invest Ophthalmol Vis Sci* 2018; 59(3):1652-1658.
 4. Sha P, et al. (Durbin M) *Invest Ophthalmol Vis Sci* 2017; 58:3510

11

Commercially Available Anterior Eye OCT

David Huang, MD, PhD www.COOL-Lab.net * Pending FDA clearance for anterior eye imaging

12




13

Yan Li, PhD Elias Pavlatos, PhD


Corneal OCT Applications

CORNEAL TOPOGRAPHY WITH OCT


14



15



16



17


Motion Detection Improves Repeatability

Diopters	No Motion Detection	Motion Detection
Anterior Mean Power	0.28	0.14
Anterior Astigmatism - Cardinal	0.92	0.28
Anterior Astigmatism - Oblique	1.16	0.24
Posterior Mean Power	0.04	0.03
Posterior Astigmatism - Cardinal	0.13	0.05
Posterior Astigmatism - Oblique	0.15	0.05


Pooled standard deviations for 20 eyes from 10 participants, 5 Repeated OCT Scans

Pavlatos E et al. *Biomed Opt Express*. 2020; 11(12):7343-7356

18



19




20

Yan Li, PhD Elias Pavlatos, PhD


Corneal OCT Applications

DIFFERENTIATING KERATOCONUS FROM WARPAGE, DRY EYE, AND OTHER IRREGULARITIES


21



22




23



24

Epithelial Concentric Thinning

- Epithelial thickness map has more than two color-scale step changes (>5 μm) inside the analytic zone
- At least one complete ring around the thinnest point




Yang Y, et al. (Li Y) J Cataract Refract Surg 2021;47(6):759-766

25

Case 1 – Normal

Axial Power Pachymetry Epithelial Thickness




Yang Y, et al. (Li Y) J Cataract Refract Surg 2021;47(6):759-766

26

Case 2 – Subclinical Keratoconus

Axial Power Pachymetry Epithelial Thickness




Yang Y, et al. (Li Y) J Cataract Refract Surg 2021;47(6):759-766

27

Case 3 – Forme Fruste Keratoconus

Axial Power Pachymetry Epithelial Thickness



Yang Y, et al. (Li Y) J Cataract Refract Surg 2021;47(6):759-766

28

Coincident Concentric Thinning Decision

Step 1
Map parameters

Pachymetric Min < 515 μm **OR**,
Pachymetric Min-Max < -71 μm **OR**,
Pachymetric SN-IT > 28 μm **OR**,
Epithelial Std Dev > 1.9 μm

YES
NO

Step 2
Map patterns

Coincident **AND**
concentric epithelial
thinning?

YES
NO

Not keratoconus

YES **Keratoconus**

NO **Not keratoconus**

Yang Y, et al. (Li Y) J Cataract Refract Surg 2021;47(6):759-766

29


Classification Accuracy of Decision Tree

- **Manifest keratoconus** – abnormal topography and CDVA < 20/20
- **Subclinical keratoconus** – abnormal topography and CDVA ≥ 20/20
- **Forme Fruste keratoconus (FFK)** – fellow eyes of asymmetric keratoconus with normal or borderline topography and CDVA ≥ 20/20


Normal (54 eyes)	Manifest Keratoconus (91 eyes)	Subclinical Keratoconus (12 eyes)	Forme Fruste Keratoconus (19 eyes)
100%	97.8%	100%	73.7%

Yang Y, et al. (Li Y) J Cataract Refract Surg 2021;47(6):759-766


30




31



32



Yan Li, PhD




Elias Pavlatos, PhD


Corneal OCT Applications

QUANTITATIVE INDICES TO SEPARATE PRIMARY AND SECONDARY EPITHELIAL CHANGES

33



34




35

Epithelial PSD is a very sensitive detector of keratoconus

	Keratoconus N= 35 eyes 150 control	20/20 Keratoconus N=50 eyes 150 control	Forme fruste keratoconus N=8 150 control
Sensitivity	100%	96.0%	87.5%
Specificity	100%	100%	100%

1. Li Y, et al. (Huang D) *Ophthalmology* 2012; 119:2425-33.
 2. Li Y, et al. (Huang D) *J Cataract Refract Surg* 2016; 42(2):284-95.

36



37

Patient Groups

Keratoconus Groups

- Manifest Keratoconus (n = 89)**
 - Topographic signs of keratoconus
 - Corrected visual acuity < 20/20
- Subclinical Keratoconus (n = 16)**
 - Topographic signs of keratoconus
 - Corrected visual acuity $\geq 20/20$
- Forme Fruste Keratoconus (n = 26)**
 - Normal topography
 - Corrected visual acuity $\geq 20/20$
 - Keratoconus in fellow eye

Non-Keratoconus Groups

- Contact Lens Warpage (n = 18)**
 - Suspicious topographic steepening
 - Reduction in corrected visual acuity
- Normal (n = 32)**
 - Normal topography
 - Corrected visual acuity $\geq 20/20$

www.COOLLab.net

38

Classification Accuracy of EM Index

Binary classification


- Class 1 = non-keratoconus (normal and warpage)
- Class 2 = keratoconus (manifest, subclinical, forme fruste)

EM Index	Classification Accuracy (%)				
Cutoff	Normal	Warpage	Manifest Keratoconus	Subclinical Keratoconus	FF Keratoconus
1.39 ± 0.01	100 ± 0	98.9 ± 2.2	100 ± 0	100 ± 0	51.5 ± 1.9

• Cutoff determined at 50% probability by logistic regression
• 5-fold cross-validation repeated 5 times

www.COOLLab.net

39



40

Concept Applies to Other Conditions

Primary Deformation


- Warpage
- Dry Eye
- EBMD

Secondary Modulation


- Keratoconus/ectasia
- Post-LASIK/PRK
- Stromal scars
- Stromal/Bowman dystrophies

www.COOLLab.net

41



Yan Li, PhD




Clara Llorens Quintana, PhD


Corneal OCT Applications

TOTAL (NET) CORNEAL POWER AND ASTIGMATISM

42



43



44

OCT Corneal Power Repeatability: Published Results

		RTVue ¹	RTVue ²
Corneal power (D)	Total (Net)	0.19	0.10
	Anterior	0.19	0.11
	Posterior	0.02	0.02

1. Tang M, et al. (Huang D) *J Cataract Refract Surg*. 2010;36(12):2115-2122.
 2. Wang Q, et al. (Huang J) *Cornea*. 2015;34(10):1266-1271.

45


Repeatability: RTVue vs Avanti

		RTVue	Avanti
Corneal power (D)	Total (Net)	0.18	0.14
	Anterior	0.20	0.15
	Posterior	0.04	0.04


n = 24 eyes of 12 normal subjects

David Huang, MD, PhD www.COOLab.net

46



47



48

Subjects

60 pseudophakic eyes (39 patients)

- With non-toric monofocal IOL
- Post Laser visual correction

↓

Manifest refraction astigmatism = ground truth

Llorens-Quintana C, et al. (Li Y) *Invest Ophthalmol Vis Sci* 2021;62:2026

49

Net Corneal Astigmatism Repeatability


Coefficient of repeatability
Pooled standard deviation* 1.96 * sqrt(2)

	Avanti*	Pentacam	Significance**
Cardinal (D)	0.22	0.50**	$p < 0.05$
Oblique (D)	0.19	0.44**	$p < 0.05$
Vector (D)	0.29	0.67**	$p < 0.05$

* Experimental software, not FDA-cleared
** F test comparison with OCT Net astigmatism

Llorens-Quintana C, et al. (Li Y) *Invest Ophthalmol Vis Sci* 2021;62:2026

50



51

OCT Corneal Mapping

- Four maps
 - Epithelial thickness
 - Pachymetry
 - Anterior topography
 - Posterior topography
- Distinguish ectasia from warpage
- May provide more accurate net corneal astigmatism measurement in aberrated corneas
- Advanced features still FDA pending

52

Grants & Material Supports

NIH Grant R01 EY029023, R01 EY028755, P30 EY010572

Unrestricted grant from Research to Prevent Blindness

Material support from Optovue, Inc.

53

www.COOLLab.net

CENTER FOR OPHTHALMIC OPTICS & LASERS

54

Measuring the cornea for IOL power calculations


Li Wang, MD, PhD,
Cullen Eye Institute
Baylor College of Medicine,
Houston, Texas, USA

Consultant:
Alcon
Carl Zeiss Meditec

1

Corneal power measurements

- ✦ Range of devices:
 - ✦ 4 points: keratometer
 - ✦ 6-32 points: ocular biometers
 - ✦ >500 points: topographic / tomographic values averaged over the central 3-4 mm zone



2

Corneal power measurements

- ✦ Most of the commonly used devices measure only the anterior corneal curvature
- ✦ Total corneal power calculated
 - ✦ Assume constant ratio of posterior to anterior corneal curvature
 - ✦ Standard index of refraction 1.3375 (USA) or 1.332 (Europe)

3

Ratio of posterior/anterior radii of curvature from Galilei*


	Ratio	"Standardized" index of refraction	Ratio in RK eyes (n=114) with Avanti OCT Mean: 0.94 Range: 0.66 – 1.30
Gullstrand's schematic eye	0.883	1.3325	
Normal (n=94)	0.816	1.3278	
Myopic-LASIK/PRK (n=61)	0.765	1.3246	
Hyperopic-LASIK/PRK (n=9)	0.857	1.3302	

* Wang L, Mahmoud AM, Anderson BL, Koch DD and Roberts CJ. IOVS 2011; 52:1716-22

4

Ratio of posterior/anterior radii of curvature

- ✦ Effect on calculation of total corneal power
 - ✦ Small in normal eyes
 - ✦ Key issue in post-refractive eyes




Posterior corneal power: -3.5 D

5


IOL power calculation in post-refractive eyes

- ✦ Traditionally using corneal power based on anterior corneal surface
- ✦ Formulas/methods using data from both anterior and posterior corneal surfaces emerging
 - ✦ OCT-based IOL calculation formula
 - ✦ Galilei TCP formula
 - ✦ Potvin-Hill Pentacam
 - ✦ Ray-tracing

6



7



8

What is ideally the best IOL formula?

- Ray tracing:
 - Incorporates all surfaces and aberrations
 - Cornea
 - IOL
 - Major limitation
 - Does not provide a solution to ELP
 - Outcomes with ray tracing formulas to date are not better than our other best formulas

9

Outcomes reported in literature: recent large studies including newer formulas

	% ± 0.5 D
Myopic LASIK/PRK eyes	40% - 85%
Hyperopic LASIK/PRK eyes	38% - 73%
RK eyes	29% - 62%

- Reported more accurate formulas:
 - Average IOL power, OCT-based, Barrett True K, Haigis-L, Maskit

10

Current best practices

- Use as many formulas as possible
- Rely more on
 - Average IOL power
 - OCT-based, Barrett True K No History, and Haigis-L
- Warn patients of inaccuracy and possible additional surgeries and costs

11

"Holy Grail": postop adjustable IOL?

- RxSight: light adjustable lens
 - Curvature change
- Perfect Lens and Clerio
 - Localized refractive index change

12

Posterior corneal astigmatism


13

Magnitude of posterior corneal astigmatism


- ◆ Our paper: 0.30 ± 0.15 (0 – 1.10) (n = 715)
- ◆ Tonn...Kohlen: 0.33 ± 0.18 (0 – 1.35) (n = 3818)

Koch et al. JCRS 2012; 38:2080-2087
Tonn B, Klaproth OK, Kohlen T. IOVS 2014; 56:291-8


14




15



16



17



18

Toric IOL calculation: how do we incorporate posterior corneal astigmatism (PCA)?

- Use regression/theoretical models
 - Baylor nomogram
 - Abulafia-Koch: Vector version of Baylor nomogram + clinical data
 - J&J (AMO): J&J clinical trial data + Baylor nomogram
 - Barrett toric calculator (standard, predicted PCA)
- Measure the posterior cornea
 - Barrett toric calculator (new, measured PCA)

19

Baylor Toric IOL Nomogram, Version 2

Temporal clear corneal incision and target for postop astigmatism of 0.4 D WTR to account for ATR shift with age

Effective IOL cylinder power at corneal plane (D)	WTR (D)	ATR (D)
0	≤ 1.69 (>1.0: PCR)	≤ 0.39
1.00	1.70 - 2.19	0.40* - 0.79
1.50	2.20 - 2.69	0.80 - 1.29
2.00	2.70 - 3.19	1.30 - 1.79
2.50	3.20 - 3.79	1.80 - 2.29
3.00	3.80 - 4.39	2.30 - 2.79
3.50	4.40 - 4.99	2.80 - 3.29
4.00	5.00 -	3.30 - 3.79

*Especially if specs have more ATR

20

Barrett toric calculator (apacrs.org)

21

Problem with regression approaches for selecting toric IOLs: Scatter

- Wide individual variation
- Even "normal" corneas

22

Does measuring the posterior cornea improve toric outcomes?


- Prior Scheimpflug studies:
 - No benefit to direct measurement
- Studies:
 - Regression models vs. direct measurements
 - More accurate with regression models

23


Our study

- To compare the accuracy of predicting residual astigmatism after cataract surgery using Barrett toric calculator with:
 - Predicted PCA
 - Measured PCA from IOLMaster 700
- Dataset from VERACITY surgical database
 - Included eyes with monofocal non-toric IOLs


24




25




26




27




28



29




30




31

- ### Toric IOL in LASIK/PRK/RK eyes
- ✦ High expectations following cataract surgery
 - ✦ Good uncorrected visual acuity
 - ✦ Spectacle independence
 - ✦ Corneal astigmatism common
 - ✦ LASIK/PRK performed to eliminate ocular refractive error including astigmatism
 - ✦ Residual or induced to compensate for lenticular astigmatism


32



33



34



35

- ### IOL toricity selection
- ✦ Posterior corneal astigmatism
 - ✦ Difficult to estimate
 - ✦ Estimate 0.3 D of against-the-rule refractive effect
 - ✦ Target for corneal astigmatism correction
 - ✦ 0.3 D WTR
 - ✦ Be conservative and not over-correct!

36

IOL spherical power selection: ascrs.org

Post-refractive IOL calculator

Version 4.3
Made possible by an unrestricted educational grant from Alcon Laboratories and The ASCRS Foundation

37

When corneas met all 3 criteria:

- % of eyes with refractive astigmatism ≤ 0.5 D
 - Post-myopic LASIK (n = 56): 80%
 - Post-hyperopic LASIK (n = 19): 84%
 - Post-RK (n = 41): 76%

Cao D, Wang L, Koch DD. JCRS . 46(4):534-539, 2020.

38

Barrett True-K toric calculator

Myopic LASIK
2D @ 90
Net astigmatism: 1.39D@92

39

Case #1

- 75 yo male
- s/p LASIK
- Cataract OS

40

Preop corneal astigmatism

MR	-0.75 x 1.0 x 175
Lenstar	1.29 @ 169
IOLMaster 700	1.30 @ 169
Galilei	
SimK	0.99 @ 159
TCP	1.37 @ 165

- 19.0 D ZCT225@172
- POM#1
- UCVA = 20/20
- MR: plano

41

Case #2: previous RK

- 60 yrs male cataract surgery OS
- MR: + 1.75 + 0.75 x 166
- 8 radial cuts and one superior T cut
- Target for D

42

Corneal astigmatism

- + MR: $+1.75 + 0.75 \times 166$
- + IOLMaster: 0.93 D @150°
- + Lenstar: 1.04 D @149°
- + Galilei SimK: 1.01 D @157°
- + Galilei TCP: 1.19 D @173°
- + Atlas: 0.76 D @177°

Toric or nontoric IOL? ZCT150 @160°

43

IOL spherical power: which would you choose??

Calculator formulas for target -0.5 D:
25.71 to 27.45 D

44

26.5D ZCT150 @160 implanted targeting -0.50

POM #1:
UCVA = 20/100
-1.75 Sph = 20/20

POM #2:
UCVA = 20/50
-1.25 Sph = 20/20

POM #3:
UCVA = 20/70
-1.75 Sph = 20/25

What would you do now?

45

Fit with -1.25 D CL

- + Realized how nice a little near vision is
- + Wants perfect distance vision OD!
- + It has 12 radial and 8 T cuts....

46

Summary

- + Accurate total corneal power estimation is crucial
 - + Posterior corneal power measurements improve accuracy, especially in un-usual eyes
- + IOL power calculation in post-refractive eyes still a ways to go
 - + Especially in RK eyes
 - + More accurate corneal power measurements and IOL power formulas are needed

47

Summary

- + Considering posterior corneal astigmatism in toric IOL selection improves accuracy
- + Toric IOLs can work well in post-refractive eyes
 - + Corneas met all 3 criteria
- + Postop IOL power adjustment is very promising

48